



OST 100/150

COMPREHENSIVE REFRACTION UNIT
FOR SMALL SPACE



OST 100/150 COMPACT SERIES

OST 150* additional features:

Chart projector holder with variable lighting



Lift-up elbow rest



Foot rest



OST 100 «STANDARD PACK» CONFIGURATION:

2 instruments sliding table
Motorized chair
Projector plate
Self-balancing refractor arm.
Led lighting

OST 150 CONFIGURATION:

Standard pack + Lift-up elbow rests
and footrest

OPTIONS:

Additional projector plate
110V power supply
Stand for keyboard
Fixed chinrest

Side unit body covers:

- Red (RAL 3001)
- Blue (RAL 5005)
- White (RAL 9016)

SPECIFICATIONS

Chair maximum load:
150 kg

Sliding tabletop:
Dimensions: 42 x 82 cm
Max. weight: 50 kg


Versions:
Right or left handed version

Power supply:
230V-50/60Hz


Dimensions:
1.25 (L) x 1.5 (P) x 1.80 (H)m

Weight:
195 kg

OST 100. Side view



OST 100. Top view



As improvements are made, images and specifications are not contractually binding and may be changed without prior notice.



ESSILOR INSTRUMENTS

81 boulevard Jean-Baptiste Oudry
94000 Créteil
France
Tél. 33(0)1 49 80 62 80
www.essilor-instruments.com

