

MEASUREMENTS MADE EASY



eye-ruler²

DIGITAL MEASUREMENTS ON IPAD



eye-ruler²

OPT FOR THE EXPERTISE OF AN ESSILOR EXCLUSIVE MEASUREMENTS SOLUTION




eye-ruler™ 2 is a compact iPad®-based measurement solution providing very accurate centration parameters for all lens types.

Benefit from the frame comparison module for ECPs, allowing customers to get inspiration and choose their favourite frame.




COMPARE AND CHOOSE YOUR CLIENT'S FAVOURITE FRAMES



FRAME CHOICE

- Select and compare up to 4 frames.
- Send the pictures to your patient by email.



ENJOY AN INTUITIVE AND EFFICIENT USER EXPERIENCE



FAR VISION MEASUREMENTS

- Pupillary distances
- Fitting heights
- Segment heights for bifocals
- Eye-lens distance
- Pantoscopic tilt
- Wrap angle
- Frame dimensions
- Effective lens diameter




NEAR VISION MEASUREMENTS

- Pupillary distances
- Fitting heights
- Reading distance

Essilor exclusive measurements:

- Near Vision Behavior parameter for Varilux® X and XR series™ lenses



MEASUREMENTS RESULTS

- User-friendly display with all results at a glance.
- 2 displays available: chart and grid.



COMPACT IN DESIGN AND USABILITY

With **eye-ruler™ 2**, opt for a comfortable user experience on a hand-held device to intuitively use the different centration parameters.

Discover the accuracy of the measurements, calculated using 3D reconstruction technology on the customer's frames position in relation to their eyes.



QUICK & EASY TO INSTALL



**WIDE COMPATIBILITY
WITH IPAD®**



ULTRA-LIGHT DEVICE



EASE OF USE



**LATEST EMBEDDED
TECHNOLOGIES (Bluetooth® low energy)**



©Essilor International – January 2024. Essilor® and eye-ruler™ 2 are trademarks of Essilor International. ©2024 Apple Inc. All rights reserved. Apple, the Apple logo, iPad® are trademarks of Apple Inc, registered in the U.S. and other countries. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth® SIG, Inc. and any use of such marks by Essilor International is under license.

ESSILOR INTERNATIONAL © - RCS Créteil 439 769 654 - M01522-01-en-Ey-ruler 2 Brochure-device-Standard - January 2024



ESSILOR INTERNATIONAL

147 rue de Paris, 94220 Charenton-le-Pont - France

Tel: +33 (0)1 49 80 62 80

www.essilor-instruments.com

